[Company name]

How to use this document: delete this paragraph before saving, sharing, or circulating it. Every part shown in blue is yours to replace: type your own details over each blue [field] until no blue remains. Add or remove table rows to fit your needs. Then select all (Ctrl+A) and set the font colour to the same dark colour throughout. What is left is your finished document.

Dental practice sterilization log

Autoclave cycles - every load logged, every indicator checked.

Date	[Date]

Sterilization log

Date	Cycle no.	Load contents	Temp / time	Indicator pass	Initials
[Date]	[no.]	[e.g. handpieces, mirrors]	[134°C / 3.5 min]	[Y/N]	

Sign-off

Name	Signature	Date